i

Krekeler Messschieber

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Krekeler caliper

Inserting Implants gets much easier and effective with the caliper acc. to Prof. Dr. med. dent. Gisbert Krekeler. Working is safer, faster and more precise by the combination of several functions in one instrument.

Functionality of the Krekeler Caliper

For a maximum stability always the widest possible implant diameter should be chosen. Due to the exactly 1.5 mm measuring rods, you can read the maximum diameter directly. Due to scientific research in order to avoid bones resorbition after the implantological treatment, there should be a wall of at least 1.5 mm of bone material to each side.

1. During the measurement of the interdental gap, it is possible to read directly the maximum implant diameter from the IN mark. A correct measurement is performed by applying the rods at the cervix dentis of the neighboring teeth. From the OUT mark you may read the the overall measure.
2. The sharpened rods also support the operator in placing the first mark in the bone, to ensure a safe distance to the neighbouring implant / tooth. If there is a neighbour implant it is possible to use one rod as spacer for the minimum distance of the next drill.
3. At the distal ending of the Krekeler Caliper you also find a drilling guide and a 5 and 7 mm graduation to measure the patients alveolar ridge from crestal. This supports an easy and direct evaluation if there is enough bone material left.