



SCHWERT
SWORD · ESPADA
COUTEAU · SPADA

SAFE-EX 2

patented



1033-61

Sicherheitsklingenhalter,
ergonomische Form,
zerlegbar
Safety blade holder,
ergonomic shape,
dismountable



**Nie wieder
Verletzungen
·
Never again
injuries**



SAFE-EX

Das alltägliche Verletzungsrisiko beim Aufnehmen und vor allem dem Entfernen der Klinge vom Skalpell Griff ist mit herkömmlichen Instrumenten unvermeidlich. Mit dem neuen SCHWERT SAFE-EX hat man eine Innovation geschaffen die den Praxisalltag verändert. Die im Griff integrierte Funktion, ähnlich der Benutzung eines Kugelschreibers, reduziert die Verletzungsgefahr auf ein Minimum. Die Klinge wird beim SAFE-EX per Knopfdruck am distalen Ende einfach ausgeworfen. Praxispersonal wird im Zuge der Wiederaufbereitung wirksam vor Krankheitsübertragungen und Infektionen durch Schnitte geschützt.

SAFE-EX

The daily risk of injuries by mounting and especially dismantling of blades is not avoidable by standard blade holders. With the new SCHWERT SAFE-EX an innovation was developed that changes the daily routine eminently. The integrated function that works like a regular ball pen reduces the risk of an injury to a minimum. The blade is ejected from the SAFE-EX just by pushing the button on the distal end. The practice staff, especially during the re-processing, is safe from infections caused by unintended cuts.



SCHWERT
SWORD · ESPADA
EPEE · SPADA

i

Einfach und unkompliziert!

Durch den neuen speziellen Schliff, ist die Klinge von beiden Seiten anbringbar.



Easy and uncomplicated!

Due to the new unique cut, the blade can be attached from both sides.



1033-51







Edelstahl, gebogen, ergonomische Form
Stainless steel, curved, ergonomic shape



1033-50

Edelstahl, gerade, ergonomische Form
Stainless steel, straight, ergonomic shape

Passende Skalpellklingen · Compatible scalpel blades

	100 Stk., steril* 100 pcs., sterile*
 FIG. 10	1063-10
 FIG. 11	1063-11
 FIG. 12	1063-12
 FIG. 12D	1063-12D
 FIG. 15	1063-15
 FIG. 15C	1063-15C