



i

SCHWERT Design · LUCAS Küretten · LUCAS curettes

Die bewährten Lucas Küretten im SCHWERT Design Griff für ein optimales Handling. Die Größen 1,5 mm, 2,5 mm und 3,0 mm sind nun neben der Standardausführung auch mit gezahnten Arbeitsenden verfügbar. Die feine Zählung ermöglicht dabei eine ideale Behandlung, im Bereich der Zystenentfernung und wirkt effektiv einem Wegrutschen, an der zu behandelnden Stelle entgegen.

The proven Lucas curettes in the SCHWERT design handle for a precise handling. The sizes 1.5 mm, 2.5 mm and 3.0 mm are now beside the standard figures also available with serrated working tips. The fine serration allows an ideal treatment in the area of removing cysts and avoids effectively a slipping away from the spot which has to be treated.



1372-0 LUCAS Kürette, Ø 1,5 mm, ungezähnt
LUCAS curette, Ø 1.5 mm, non-serrated



1372-1 LUCAS Kürette, Ø 2,5 mm, ungezähnt
LUCAS curette, Ø 2.5 mm, non-serrated



1372-2 LUCAS Kürette, Ø 3,0 mm, ungezähnt
LUCAS curette, Ø 3.0 mm, non-serrated



1372-10 LUCAS Kürette, Ø 1,5 mm, gezähnt
LUCAS curette, Ø 1.5 mm, serrated



1372-11 LUCAS Kürette, Ø 2,5 mm, gezähnt
LUCAS curette, Ø 2.5 mm, serrated



1372-12 LUCAS Kürette, Ø 3,0 mm, gezähnt
LUCAS curette, Ø 3.0 mm, serrated